PROVA N. 1

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
https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf:

Fume Hood Safety. Test your fume hood containment before first use and at least once per year thereafter. Verify the exhaust system is operating before each use.

Handwritten signature and initials in black ink, appearing to be 'M. R. M. D. 9/6'.

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(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

Verify operation of the fume hood alarm before each use. This may be done by temporarily blocking the air sampling port to simulate loss of airflow. Close the sash as much as possible when working at the fume hood - use the glass as a face shield.

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
(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

Hoods equipped with a combination vertical and horizontal sash are intended to be used primarily with the vertical frame down and the horizontal moving panels open. Never raise the sash vertically without closing the horizontal panels first.

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Close the sash completely when not working at the fume hood. Avoid placing your head inside the fume hood. Avoid rapid movements in front of and inside the fume hood. This can cause air current to draw contaminants out of the fume hood.

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Avoid raising the sash(s) too quickly. Raise the sash at a moderate rate and allow the air velocity to increase or the variable air volume system to react to the door movement. Never work with sashes set higher than the established working height. Contact safety or air balancing staff if the safe working height is unknown.

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Never work with both sashes open on a walk-in hood. The lower sash is for setup access only. Turn on the lights. Good illumination will reduce eyestrain. Set up apparatus as close to the back as possible. In no case should apparatus be closer than 6 inches (15 cm) from the door.

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Don't block airflow. Raise large objects 2 inches (5 cm) off the counter by placing them on blocks. This allows airflow underneath and prevents stagnant areas. Never store large quantities of chemicals in the fume hood. Chemical containers block airflow and create unnecessary hazards. Use cabinets designed to store chemicals instead.

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(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

Clean up chemical spills as soon as possible. Avoid damage to countertops and buildup of large quantities of vapors. All electrical equipment connected to fume hood receptacles must be certified or approved and in good condition including electrical ground if applicable.

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(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

Avoid excessive heat within the fume hood. Special hood configurations are available for high heat situations.

Baffles should be adjusted by qualified air flow experts or safety officers only.

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(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

Perchloric acid hoods should have the wash down system operated at least once per day until a thorough cleaning is accomplished. 5-15 min will probably be sufficient depending upon flow rates achieved by the in-duct portion of the system.

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
(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

Never allow the area immediately in front of a fume hood in use to become a traffic area. Others walking by the hood face will adversely affect containment. When inserting and removing arms and objects into the fume hood, ensure that movement is parallel with the airflow direction. Avoid any scooping type motions. Avoid rapid removal of objects or arms from hood.

A handwritten signature in black ink, appearing to be 'MR. M. D. Y.', located at the bottom left of the page.

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
Avoid filling the hood with excessive equipment or materials. As a rule, no more than 50% of the available work surface should be covered with materials or equipment.

A handwritten signature in black ink, appearing to be 'M. C. M. Z. 1/16', located at the bottom left of the page.

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
Fume Hood Alarms. Your fume hood must be equipped with a suitable low flow alarm*. The low flow alarm is designed to protect the user from exposing themselves to hazards by alarming when low face velocity is detected. Do not mute the alarm or otherwise defeat its operation. Such alarms may be factory supplied or purchased separately.

Handwritten signature and initials, possibly 'MR' followed by a name and a date '10/16'.

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
(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

When the alarm sounds immediately close the sash and step away from the hood. If the alarm does not reset itself, evacuate the area and contact facilities. There may be a malfunction in the exhaust system serving that hood.

The bottom of the page contains several handwritten signatures and initials in black ink. On the left, there is a large, stylized signature that appears to be 'MR'. To its right is another signature, followed by a set of initials 'JR' and a date '4/16'.

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Under certain circumstances, it may be necessary to raise the sash above the established working height to allow materials and apparatus to enter the hood. In this case, wait two minutes after all hazard generating activity has ceased and all containers of materials are tightly closed.

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
(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

The sash may then be raised as high as needed for setup operations. During this time only, it is safe to temporarily mute the alarm by pressing the enter button.

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
(https://www.newenglandlab.com/userfiles/files/resources/Product%20Manuals/Fume_Hoods/Fume_Hood_User_Installation_Product_Manual.pdf):

Calibration: Generally, low flow alarms must be calibrated in the field after the fume hood exhaust system has been adjusted or balanced. This service is not included with the purchase of the hood and it is recommended that this be performed by a firm specializing in air balancing or ventilation systems.

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
Written procedures should be developed by your safety officer or committee to handle a variety of emergencies such as: Fire, Explosions, Out of control experiments, Power Failures, Blower Failure

The bottom of the page contains four handwritten signatures or initials in black ink. From left to right, they appear to be: a stylized 'PK', a signature that looks like 'Anna', a signature that looks like 'D', and a signature that looks like 'YB'.