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This rheometer enables the user to perform rheological measurements in CR (Controlled Rate), CS (Controlled Stress), and CD (controlled Deformation) mode in rotation and in CS and CD modes in oscillation.

In the case of a rotational rheometer, the viscosity of a liquid is calculated in accordance with the Newtonian conditional equation for viscosity:

Viscosity (h) = Shear Stress (t) / Shear Rate (g)

at defined ambient conditions regarding measuring time, temperature and pressure. In rheometers, operating in accordance with the CR-principle, a speed (angular velocity) is preset which, in the sensor system filled with a sample, causes a shear rate. The torque required for achieving and maintaining the desired shear rate is the viscosity-proportional parameter.

CS-rheometers are designed to operate according to the reversed principle. Here a torque (shear stress) is present and the resulting movement (deformation), i.e. the resulting angular velocity (shear rate) is measured. The measurement with rotational rheometers can be summed up to the predetermination of a force and from the measurement of the resulting movement, a suitable geometry can be derived.

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Time: Have a quick look at some intense peaks and adjust the time to get a decent signal/noise ratio.

Increment: It is recommended to choose the 2Theta step such that you obtain at least 7-8 points above the full width at half maximum of the Bragg peak in order to get good peak shapes. This is especially important if you want to model data yourself or have them modelled by someone else (profile fitting, Rietveld refinement...). Anything above 12 points is usually a waste of measurement time.

PSD opening: (Detector slit) Can be at the maximum. As with the slit, play around for resolution.

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The spectrum measurement is mainly used in qualitative analysis of the sample. Its direct diagram display mode makes the character of the sample clear and provides you the absorption degree of the light of the material at certain wavelength point or within certain wavelength range, making it an important function of the UV-VIS spectrophotometer.

Quantitation Measurement: In quantitation measurement, the measurement value of the sample to be measured is compared to the standard curve to calculate the concentration of the sample to be measured, namely draw a absorbance-concentration relation curve according to the absorbance of the sample with known concentration, then measure the Abs value of the sample to be measured under the same conditions and calculate the concentration of the sample to be measured according to the Lambert-Beer law (the absorption of the material of the light is positively proportional to the concentration of the material) and the standard curve.

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